




# INTREPED<sup>®</sup>

Intraosseous Fusion Device



## Advancing Fixation

New era of TMT fusion technology.

## Streamlined Procedure

Less hardware, simplified joint preparation.

## Smaller Incision

Compared to conventional TMT fusion procedures.

## No Dorsal Prominence

Lower profile than traditional plating methods.

*Step into the Future with INTREPED<sup>®</sup>*



**AUXANO MEDICAL**





# INTREPED®

## Intraosseous Fusion Device

### INTREPED® – A Step Forward in Fusion

The INTREPED® addresses the painful arthritic midfoot in a novel and efficient manner. The unique, interrupted thread geometry provides initial fixation stability, while the fenestrations enable bone growth throughout the titanium implant. Patient healing & comfort are promoted with smaller incisions and recessed implant position.



### INTREPED® System

01-0001 – KIT CONTENTS									
IMPLANT – Qty 2 ea		DRILL BIT – Qty 2 ea		DRIVER			SPANNER DRIVE		
SIZE (mm)	ORDER #	SIZE (mm)	ORDER #	SIZE (mm)	QTY	ORDER #	SIZE (mm)	QTY	ORDER #
06 x 10	01-10610	06	01-40006	06	2	01-40206	N/A	1	01-50003
06 x 15	01-10615	09	01-40009	09	2	01-40209	K-WIRE		
06 x 20	01-10620						1.6 x 6	12	01-50160
09 x 10	01-10910	12	01-40012	12	1	01-40212	T-HANDLE		
09 x 15	01-10915	TAP – Qty 1 ea		SIZER			T-HANDLE		
09 x 20	01-10920	06	01-40106	06	1	01-40306	N/A	1	01-50002
12 x 10	01-11210	09	01-40109	09	1	01-40309	AXIAL HANDLE		
12 x 15	01-11215						N/A	1	01-50001
12 x 20	01-11220	12	01-40112	12	1	01-40312	N/A	1	01-50001

The INTREPED® system is supplied non-sterile, and comes in a steam-sterilizable tray.

For more information visit us at [www.AuxanoMedical.net](http://www.AuxanoMedical.net)  
**440 262 2000 | CS@AUXANOMEDICAL.NET**



AuxanoMedical.net



Auxano Medical  
 8006 Katherine Boulevard  
 Brecksville, OH 44141  
[www.AuxanoMedical.net](http://www.AuxanoMedical.net)  
 Made In USA



© 2024 Auxano Medical 01-60014-C

\*This document is not a substitute for the IFU/DFU (Instructions for Use/Directions for Use).