


INTREPED[®]

Intraosseous Fusion Device



Advancing Fixation

New era of TMT fusion technology.

Streamlined Procedure


Less hardware, simplified joint preparation.

Smaller Incision

Compared to conventional TMT fusion procedures.

No Dorsal Prominence

Lower profile than traditional plating methods.



Step into the Future with INTREPED[®]



AUXANO MEDICAL



INTREPED®

Intraosseous Fusion Device

INTREPED® – A Step Forward in Fusion

The INTREPED® addresses the painful arthritic midfoot in a novel and efficient manner. The unique, interrupted thread geometry provides initial fixation stability, while the fenestrations enable bone growth throughout the titanium implant. Patient healing & comfort are promoted with smaller incisions and recessed implant position.


INTREPED® Technology

Central Graft Cavity

Threaded Geometry

Thru-Growth Fenestrations

AUXANO® In-Growth Technology



INTREPED® Intraosseous Fusion Device

IMPLANT SIZES – PART NUMBERS					
Diameter	Length				
	10 mm	12 mm	14 mm	15 mm	20 mm
6 mm	01-10610	–	–	01-10615	01-10620
8 mm	01-10810	01-10812	01-10814	–	–
9 mm	01-10910	01-10912	01-10914	01-10915	01-10920
10 mm	01-11010	01-11012	01-11014	–	–
12 mm	01-11210	–	–	01-11215	01-11220

The INTREPED® system is supplied non-sterile, and comes in a steam-sterilizable tray. It contains the following instruments: k-wires, depth gauge, sizers, drills, drill stops, taps, drivers, spanner drive and handle.

For more information visit us at www.AuxanoMedical.net
440 262 2000 | CS@AUXANOMEDICAL.NET



AuxanoMedical.net



Auxano Medical
8006 Katherine Boulevard
Brecksville, OH 44141
www.AuxanoMedical.net
Made In USA

