

INTREPED[®]

Intraosseous Fusion Device



Advanced Fixation

Optimized for rapid, direct integration and early column bridging.

Simple Procedure


Minimally invasive, soft-tissue sparing, faster recovery.

Zero Prominence

Avoids soft-tissue irritation common with plates and staples.

Anatomy Sparing

Simple joint preparation that preserves column length.



AUXANO[®] MEDICAL
Step into the Future[™]





INTREPED[®]

Intraosseous Fusion Device

INTREPED[®] – A Step Forward in Fusion

The INTREPED[®] addresses the painful arthritic midfoot in a novel and efficient manner. The unique, interrupted thread geometry provides initial fixation stability, while the fenestrations enable bone growth throughout the titanium implant. Patient healing & comfort are promoted with smaller incisions and recessed implant position.

INTREPED[®] Technology



INTREPED[®] Intraosseous Fusion Device

IMPLANT SIZES – PART NUMBERS					
Diameter	Length				
	10 mm	12 mm	14 mm	15 mm	20 mm
6 mm	01-10610	-	-	01-10615	01-10620
8 mm	01-10810	01-10812	01-10814	-	-
9 mm	01-10910	01-10912	01-10914	01-10915	01-10920
10 mm	01-11010	01-11012	01-11014	-	-
12 mm	01-11210	-	-	01-11215	01-11220

The INTREPED[®] system is supplied non-sterile, and comes in a steam-sterilizable tray. It contains the following instruments: k-wires, depth gauge, sizers, drills, drill stops, taps, drivers, spanner drive and handle.

For more information visit us at www.AuxanoMedical.net
 440 262 2000 | CS@AUXANOMEDICAL.NET



Auxano[®] Medical
 8006 Katherine Boulevard
 Brecksville, OH 44141
www.AuxanoMedical.net
 Made In USA

